



Resuscitator KRE-A100

www.nanomed.cc info@nanomed.cc

Resuscitator KRE-A100

Resuscitator KRE-A100 is an auto-clavable, hand held, emergency device intended for pulmonary resuscitation. Inbuilt pressure relief valve minimizes the risk of over pressure. Transparent bag for complete internal visibility. Made of durable, firm and flexible material. Substantial range of reservoir volume. Easily detachable for cleaning and disinfecting.

Features:

- 100% Latex free
- Textured grip (non-slip)
- Autoclavable repeatedly at 134 °C
- Available with mask, patient valve, intake valve, resuscitation bag
- Available in adult
- 2000 ml bag reservoir volume

Applications:

Used in emergency medical services, crash carts and critical care to provide positive pressure to inflate the lungs of an unconscious person who is not breathing, in order to keep them oxygenated.

Specifications :

Bag reservoir volume	2000 ml
Dead space (patient valve)	7 ml
Dimensions (L × D)	325 × 132 mm
Patient weight	? 33 kg
Pressure	60 cmH ₂ O
Resuscitator volume	1630 ml
Stroke volume	800 ml
Type	Adult



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc