



Resuscitator KRE-A101

www.nanomed.cc info@nanomed.cc

Resuscitator KRE-A101

Resuscitator KRE-A101 is an auto-clavable, hand held, emergency device intended for pulmonary resuscitation. Inbuilt pressure relief valve minimizes the risk of over pressure. Transparent bag for complete internal visibility. Made of durable, firm and flexible material. Substantial range of reservoir volume. Easily detachable for cleaning and disinfecting.

Features:

- 100% Latex free
- Textured grip (non-slip)
- Autoclavable repeatedly at 134 °C
- Available with mask, patient valve, intake valve, resuscitation bag
- Available in pediatric
- 1600 ml bag reservoir volume

Applications:

Used in emergency medical services, crash carts and critical care to provide positive pressure to inflate the lungs of an unconscious person who is not breathing, in order to keep them oxygenated.

Specifications :

Bag reservoir volume	1600 ml
Dead space (patient valve)	7 ml
Dimensions (L × D)	262 × 105 mm
Patient weight	10 kg ~ 33 kg
Pressure	40 cmH ₂ O
Resuscitator volume	680 ml
Stroke volume	350 ml
Type	Pediatric



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc