



Electric Breast Pump KBP-A100

www.nanomed.cc info@nanomed.cc

Electric Breast Pump KBP-A100

Electric Breast Pump KBP-A100 is a dual expression high pressure electric breast pump with adjustable suction control. Suction pump comprises of feeder system integrated with suction assembly which can be sterilized. Soft suction cup design simulates natural feeding action for gentle expression. Closed one-way pumping mechanism prevents backflow or spillage.

Features:

Double-Expression technology
Latex-free Clean assembly
Comfortable expression
Soft petal massage breast shield

Applications:

Used in nursing homes, gynecology wards, NICU as well as home-care for either working mothers or to stimulate lactation in women or to relieve engorgement or aid in recovery from pregnancy.

Specifications :

Bottle : Temperature indication	Milk Temp 40 °C – yellow Milk Temp > 40 °C – Fluorescent yellow
Construction	Food grade PP and Silicone, BPA free.
Dimension	290 x 110 x 270 mm
Display	LCD display
Expression	Double-phase
Feeding Bottle	Night light with Temperature sensor Technology
Modes	Massage and Expression
Operational Levels	9 grades of adjustment
Power Supply	AC power or Battery
Pumping System	Back flow protection mechanism
Weight	1 Kg



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc