



Nebulizer KNM-A100

www.nanomed.cc info@nanomed.cc

Nebulizer KNM-A100

Nebulizer KNM-A100 is a portable, compact tabletop ultrasonic nebulizer unit. Equipped with mesh membrane nebulization technology for consistent aerosol particle formation. Control system features high nebulization flow-capacity, deeper and penetration at a short inhalation time as per safe standards of performance.

Features:

- Ultrasonic vibration technology
- High-velocity flow
- Plug and play operation
- Overflow protection
- Battery powered
- 2 Nebulizer cup sizes

Applications:

Used in hospital wards or home-care for deep penetration in respiratory system with aerosol medicinal particles for aerosol treatments reducing side effects due to lower medicine consumption.

Specifications :

Continuous Performance	0 - 4 hrs
Dimensions	290 x 200 x 250 mm
Filter Particle Dimension	4.2 microns
Max Medication Capacity	350 ml (Large cup) 150 ml (Small cup)
Nebulization Rate	4 ml/min
Power Supply	230/120 V, 50/60Hz
Timing Range	0-60 min (no level for adjustment)
Weight	2 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc