



## Electric Suction Machine KSM-A100

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

## Electric Suction Machine KSM-A100

Electric Suction Machine KSM-A100 is a light-weight, portable phlegm suction unit. Incorporated oil-free suction pump with analog vacuum gauge housed in a crevice-free plastic exterior ensures precise infection-control, fast suction build up and improves flux-rate for efficient suction performance. Advanced overflow-protection design of graduated suction bottle prevents leakage of aspirated fluids.

### Features:

- Oil-free suction mechanism
- Advanced overflow-protection design
- Compact design
- Suitable for emergency use

### Applications:

Used in hospital wards, surgical rooms as well as at home by medical professionals or patients for performing aspirations or suction of congestive secretions for respiratory home-care.

## Specifications :

|                             |                        |
|-----------------------------|------------------------|
| <b>Continuous Operation</b> | Up to 30 min           |
| <b>Dimensions</b>           | 350 x 260 x 290 mm     |
| <b>Input Power</b>          | 90 VA                  |
| <b>Noise</b>                | ? 65 dB                |
| <b>Power Supply</b>         | 220 V±10% , 50 Hz±2%   |
| <b>Pumping Rate</b>         | 15 L/min               |
| <b>Reservoir Capacity</b>   | 1000 ml (1 Jar)        |
| <b>Vacuum Pressure</b>      | ( Negative ) 0.075 MPa |
| <b>Weight</b>               | 5 Kgs                  |



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)