



Vein Illuminator KVI-A100

www.nanomed.cc info@nanomed.cc

Vein Illuminator KVI-A100

Vein Illuminator KVI-A100 is a portable illumination unit with NIR for projection of vein image on superficial skin. It allows quick and precise vein display which is useful for infants (newborns), pediatricians and adults. Provides significant help for patients whose vein access becomes a challenge such as dark skin, hairy skin, swollen tissues, obese, old age etc.

Features:

- Large display screen for vein viewing
- Provides real time accurate images
- Adjustable height antennae
- Improves site selection (cursor locator)
- Reduces discomfort

Applications:

Used in hospitals, clinical labs, blood donation camps, rehabilitation department, academic institutions and others for intravenous access and blood collection

Specifications :

Dimension	242 x 182 x 70 mm
Display Mode	Screen
Image Resolution	640 x 480 Pixel
Optimal Imaging Distance	210 mm ± 30 mm
Power Source	3000 mAh
Source Voltage	DC 7.4 V
Wavelength (Infrared)	850 nm
Weight	950 gms



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc