



Hospital bed KHB-A100

www.nanomed.cc info@nanomed.cc

Hospital bed KHB-A100

Hospital bed KHB-A100 is a standard hospital bed equipped with aluminium folding side-rails for weak or geriatric patients. Powder coated steel tube frame based with stable stationary bed base at fixed bed height renders high strength and impact resistance along combined with skid-proof rubber tip pads. ABS moulded head and foot board ensures superior strength and increased durability.

Features:

- 1 function manual bed inclinations (Fowler variations)
- Adopts crank tension mechanism
- Over-loading protection system
- Steel tube frame with S.S. head and foot board

Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

Specifications :

Construction	Electrostatic powder coated S.S. tube frame and ABS head and foot board
Dimensions	2120 x 970 x 530 mm
Load capacity	250 Kgs
Other components	Urine hooks, I.V drip, heavy duty castors
Weight	57 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc