



Hospital bed KHB-A207

www.nanomed.cc info@nanomed.cc

Hospital bed KHB-A207

Hospital bed KHB-A207 is a standard mechanical ward bed with Aluminium side-rails, overload protection hand crank and manual crank spring tensioning system ensures improved comfort and protection at various angular bed configurations for knee-rest and back-rest positions.

Features:

- 2 functional bed inclinations (Fowler/PLR variations)
- Adopts crank tension mechanism
- Over-loading protection system
- Steel tube frame with ABS head and foot board

Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

Specifications :

Back-rest angle	0° - 80° (±5°)
Construction	Electrostatic powder coated S.S. tube frame with ABS head and foot board
Dimensions	2120 x 970 x 530 mm
Knee-rest angle	0° - 40° (±5°)
Load capacity	250 Kgs
Other components	Urine hooks, I.V drip, heavy duty castors
Weight	79 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc