



## Hospital bed KHB-A211

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

## Hospital bed KHB-A211

Hospital bed KHB-A211 is a standard mechanical ward bed equipped with stainless steel guard rails and double protection screw system which improves patient safety for weak or geriatric patients along with robust stainless steel bed board which ensures superior strength and increased durability. Bed planes inclinations ensures improved comfort and protection at various angular bed configurations for knee-rest and back-rest positions.

### Features:

2 functional bed inclinations (Fowler/PLR variations)

Adopts crank tension mechanism

S.S. guard rails for support

Steel tube frame with S.S. head and foot board

### Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

### Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

## Specifications :

<b>Back-rest angle</b>	0° - 80° ( $\pm 5^\circ$ )
<b>Construction</b>	Electrostatic powder coated S.S. tube frame and head and foot board
<b>Dimensions</b>	2000 x 900 x 530 mm
<b>Knee-rest angle</b>	0° - 40° ( $\pm 5^\circ$ )
<b>Load capacity</b>	250 Kgs
<b>Other components</b>	Urine hooks, I.V drip, heavy duty castors
<b>Weight</b>	82 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)