



Hospital bed KHB-A212

www.nanomed.cc info@nanomed.cc

Hospital bed KHB-A212

Hospital bed KHB-A212 is a standard mechanical ward bed equipped with stainless steel guard rails which improves patient safety for weak or geriatric patients along with robust stainless steel fence type bed board and arc shaped head and foot board which ensures superior strength and increased durability. Bed planes inclinations ensures improved comfort and protection at various angular bed configurations for knee-rest and back-rest positions.

Features:

- 2 functional bed inclinations (Fowler/PLR variations)
- Adopts crank tension mechanism
- S.S. guard rails for support
- Steel tube frame with S.S. head and foot board

Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

Specifications :

Back-rest angle	0° - 80° ($\pm 5^\circ$)
Construction	Powder coated S.S. tube frame and arc shaped head and foot board
Dimensions	2000 x 900 x 530 mm
Knee-rest angle	0° - 40° ($\pm 5^\circ$)
Load capacity	250 Kgs
Other components	Urine hooks, I.V drip, heavy duty castors
Weight	82 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc