



Hospital bed KHB-A216

www.nanomed.cc info@nanomed.cc

Hospital bed KHB-A216

Hospital bed KHB-A216 is a standard mechanical ward bed equipped with aluminium folding side-rails, overload protection hand crank and manual crank spring tensioning system which ensures improved comfort and protection at various angular bed configurations for back-rest positions. S.S. head and foot board ensures superior strength and increased durability.

Features:

- 1 functional bed inclinations (Fowler variations)
- Adopts crank tension mechanism
- Over-loading protection system
- Steel tube frame with ABS head and foot board

Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

Specifications :

Back-rest angle	0° - 80° (±5°)
Construction	Electro-static powder coated S.S. tube frame and head and foot board
Dimensions	2120 x 970 x 530 mm
Load capacity	250 Kgs
Other components	Urine hooks, I.V drip, heavy duty castors
Weight	57 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc