



Hospital bed KHB-A218

www.nanomed.cc info@nanomed.cc

Hospital bed KHB-A218

Hospital bed KHB-A218 is a standard mechanical ward bed equipped with anti-skid moisture-proof foot cover, overload protection hand crank and manual crank spring tensioning system which ensures improved comfort and protection at various angular bed configurations for back-rest positions. S.S. head and foot board ensures superior strength and increased durability.

Features:

- 1 functional bed inclinations (Fowler variations)
- Adopts crank tension mechanism
- Over-loading protection system
- Steel tube frame with S.S. head and foot board

Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

Specifications :

Back-rest angle	0° - 80° (±5°)
Construction	S.S. tube bed board frame and head and foot board
Dimensions	2000 x 900 x 530 mm
Load capacity	250 Kgs
Other components	Urine hooks, I.V drip, heavy duty castors
Weight	60 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou,
Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc