



## Hospital bed KHB-A303

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

# Hospital bed KHB-A303

Hospital bed KHB-A303 is an electrically controlled collapsible standard ward bed with steel with powder coated side-rails for convenient patient access. Bed can be adjusted at any bed height and various angular bed configurations for back-rest and knee-rest positions. Equipped with advanced medical motors and centralized castor brake system for smooth transition at minimal effort.

## Features:

- 3 functional bed inclinations (Fowler/PLR variations)
- Advanced medical motor system
- Collapsible bed frame
- Steel frame with ABS head and foot board

## Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

## Standard accessories:

| Sr. | Standard       | Qty |
|-----|----------------|-----|
| 1   | I.V. poles     | 1   |
| 2   | Urine bag hook | 1   |

## Specifications :

|                             |                                                                          |
|-----------------------------|--------------------------------------------------------------------------|
| <b>Back-rest angle</b>      | 0° - 80° (±5°)                                                           |
| <b>Construction</b>         | Electrostatic powder coated S.S. tube frame with ABS head and foot board |
| <b>Dimensions</b>           | 2120 x 970 x 450 mm                                                      |
| <b>Height adjustability</b> | 450 mm - 700 mm                                                          |
| <b>Knee-rest angle</b>      | 0° - 40° (±5°)                                                           |
| <b>Load capacity</b>        | 250 Kgs                                                                  |
| <b>Other components</b>     | Urine hooks, I.V drip, heavy duty castors                                |
| <b>System functions</b>     | Centrally-controlled brakes                                              |
| <b>Weight</b>               | 105 Kgs                                                                  |



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)