



## Hospital bed KHB-A305

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

## Hospital bed KHB-A305

Hospital bed KHB-A305 is an electrically controlled collapsible standard ward bed with folding aluminium side-rails for convenient patient access. Bed can be adjusted at various bed angles for back-rest positions. Equipped with advanced medical motors for smooth transition at minimal effort.

### Features:

- 1-functional bed inclinations (Fowler variations)
- Advanced medical motor system
- Steel frame with ABS head and foot board
- Aluminium folding side-rails

### Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

### Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

## Specifications :

<b>Back-rest angle</b>	0° - 80° ( $\pm 5^\circ$ )
<b>Construction</b>	Electrostatic powder coated S.S. tube frame with ABS head and foot board
<b>Dimensions</b>	2120 x 970 x 530 mm
<b>Load capacity</b>	250 Kgs
<b>Other components</b>	Urine hooks, I.V drip, heavy duty castors
<b>Weight</b>	72 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)