



## Hospital bed KHB-A307

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

# Hospital bed KHB-A307

Hospital bed KHB-A307 is an electrically controlled standard ward bed with aluminium guard-rails for convenient patient access. Bed can be adjusted at various angular bed configurations for back-rest and foot-rest positions. Equipped with advanced medical motors, control panel and stainless steel portal traction frame for smooth transition at minimal effort.

## Features:

- 3 functional bed inclinations (Fowler/PLR variations)
- Advanced medical motor system
- Steel frame with ABS head and foot board
- Aluminium guard-rails

## Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

## Standard accessories:

| Sr. | Standard       | Qty |
|-----|----------------|-----|
| 1   | I.V. poles     | 1   |
| 2   | Urine bag hook | 1   |

## Specifications :

|                             |                                                                          |
|-----------------------------|--------------------------------------------------------------------------|
| <b>Back-rest angle</b>      | 0° - 80°                                                                 |
| <b>Construction</b>         | Electrostatic powder coated S.S. tube frame with ABS head and foot board |
| <b>Dimensions</b>           | 2120 x 900 x 530 mm                                                      |
| <b>Height adjustability</b> | 530 mm – 2000 mm                                                         |
| <b>Left Foot Raise</b>      | 0° - 40°                                                                 |
| <b>Load capacity</b>        | 250 Kgs                                                                  |
| <b>Other components</b>     | Urine hooks, I.V drip, heavy duty castors                                |
| <b>Right foot raise</b>     | 0° - 40°                                                                 |
| <b>Weight</b>               | 141 Kgs                                                                  |



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)