



Hospital bed KHB-A400

www.nanomed.cc info@nanomed.cc

Hospital bed KHB-A400

Hospital bed KHB-A400 is a multi-functional electrically controlled ICU bed with drop-down side-panels along with touch keypads for convenient patient access. Specialized motors enable patient positioning at any bed height and various angular bed configurations for back-rest and knee-rest positions. Equipped with one-touch intelligent bed configuration system with 10 preset positions for uncompromised patient care.

Features:

- Electrically controlled with UPS
- Multiple functional bed inclinations
- Steel frame with ABS head and foot board
- Touch-key panel control
- Data trends and storage
- ‘Patient Out-of-Bed’ alarm

Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

Specifications :

Back-rest angle	0° - 65° ($\pm 5^\circ$)
Construction	Powder-coated S.S. tube frame with ABS head and foot board
Dimensions	2330 x 1050 x 710 mm
Height adjustability	530 mm - 710 mm
Knee-rest angle	0° - 43° ($\pm 5^\circ$)
Load capacity	250 Kgs
Other components	Urine hooks, I.V drip, heavy duty castors
Reverse Trendelenburg angle	0° - 10°
System functions	Data storage, Alarm, Reset, Weighing
Trendelenburg angle	0° - 10°
User interface	Touch-key panel
Weight	148 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc