



Hospital bed KHB-A401

www.nanomed.cc info@nanomed.cc

Hospital bed KHB-A401

Hospital bed KHB-A401 is an electrically controlled ICU bed with streamline design, drop-down side-panels along with touch keypad for convenient patient access. Specialized motors enable patient positioning at any bed height and various angular bed configurations for back-rest and knee-rest positions. Equipped with one-touch intelligent patient positioning system offering 10 pre-set bed configurations for uncompromised patient care.

Features:

- Multiple functional bed inclinations
- Steel frame with ABS head and foot board
- Touch-key panel control
- Electric lock and unlock system

Applications:

Used in hospital wards, primary healthcare facilities and nursing homes with specialized functions dedicated for specific wards such as general ward, nursing, orthopedic, burns, ICU, etc. designed to be adjustable at various inclinations, heights for increased comfort and well-being of the patient.

Standard accessories:

Sr.	Standard	Qty
1	I.V. poles	1
2	Urine bag hook	1

Specifications :

Back-rest angle	0° - 65° ($\pm 5^\circ$)
Construction	Powder-coated S.S. tube frame with ABS head and foot board
Dimensions	2330 x 1050 x 700 mm
Height adjustability	450 mm - 700 mm
Knee-rest angle	0° - 43° ($\pm 5^\circ$)
Load capacity	250 Kgs
Other components	Urine hooks, I.V drip, heavy duty castors
Reverse Trendelenburg angle	0° - 10°
System functions	Reset, Electric lock and unlock
Trendelenburg angle	0° - 10°
User interface	Touch-key panel
Weight	148 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc