

Medical Emergency Trolley KTR-A100

Medical Emergency Trolley KTR-A100

www.nanomed.cc info@nanomed.cc

Medical Emergency Trolley KTR-A100

Medical Emergency Trolley KTR-A100 is a multi-purpose transport cart built with high strength steel. Robust exterior and ergonomic design features superior durability and enhanced utility. Adjustable triple-sectioned aluminium drawers provide spacious and hygienic storage solution.

Features:

- Adjustable drawer partitions
- Integrated aluminium-wire basket
- Multi-functional hooks and I.V. pole
- Integrated waste bin

Specifications :

Casters Wheels	Diagonal brake anti-static caster wheels
Construction	High strength steel
Dimensions	670 x 450 x 1000 mm
Disposal Bins	Angular-tilt ABS dispensing boxes
Drawers	Aluminium slide-out drawers
Interior Storage	3 sections with adjustable partitioning
Number of Drawers	5 standard



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou,
Guangdong, China
Email: info@nanomed.cc | Website: www.nanomed.cc