



Haemocytometer KHT-A100

www.nanomed.cc info@nanomed.cc

Haemocytometer KHT-A100

Abimed Haemocytometer KHT-A100 has an attractive appearance, shows clear and rational data construction and is convenient to use.

Features:

The Haemocytometer can be controlled through microcomputer
Displays 10 groups of numbers and can total accurately
It is easy and convenient to use

Applications:

Used in medium and small sized laboratories for counting of blood cells.

Specifications :

Count Capacity of every group	0 - 99
External Dimensions	220 x 166 x 78 mm
Maximum Count Category	10 groups
Power Rating	5 W
Voltage	220 V ± 10 %(50 Hz)



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou,
Guangdong, China
Email: info@nanomed.cc | Website: www.nanomed.cc