



High-Speed Dental Hand piece KDH-A101

www.nanomed.cc info@nanomed.cc

High-Speed Dental Hand piece KDH-A101

High-Speed Dental Hand piece KDH-A101 is a contra-angled air-driven drill instrument. Incorporated high precision gears with elongated nozzle allows convenient movements and strong torque speed ensures smooth cut. Single-point water spray jet induces on-spot cooling effect.

Features:

Seamless Scratch-resistant surface
Strong torque speed
Enhanced grip control

Applications:

Used in dental clinics for preparation of tooth for various procedures primarily suitable for cavity preparation, removal of tooth structure, margin creation and finishing restorations.

Specifications :

Air Pressure	0.2 MPa ~ 0.24 MPa
Applicable Bur Size	Diameter : 1.59 mm ~ 1.6 mm Length : 21 mm ~ 23 mm
Bur Insertion	Standard Push button
Cartridge	Collect-type
Construction	Stainless steel exterior
Drive	Air driven
Gear Ratio	1:1
Noise	68 dB
Rotation Speed	Up to 3,60,000 rpm
Spray	Single-point water spray
Torque Speed	3,00,000 rpm
Tubing Connection	2 or 4 hole
Type of Bearing	Ceramic
Weight	Approx. 180 g



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc