



## High-Speed Dental Hand piece KDH-A200

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

## High-Speed Dental Hand piece KDH-A200

High-Speed Dental Hand piece KDH-A200 is a contra-angled fiber optic dental drill instrument. Incorporated high precision spin-motor allows stable rotation ensuring high-end performance. LED light illumination and triple-point water spray improves scope of visibility.

### Features:

- Optical LED illumination
- Anti-retraction valve mechanism
- Enhanced grip control
- Quick coupling

### Applications:

Used in dental clinics for preparation of tooth for various procedures primarily suitable for cavity preparation, removal of tooth structure, margin creation and finishing restorations.

## Specifications :

|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Air Pressure</b>        | 0.2 MPa ~ 0.24 MPa                                    |
| <b>Applicable Bur Size</b> | Diameter : 1.59 mm ~ 1.6 mm<br>Length : 21 mm ~ 23 mm |
| <b>Bur Insertion</b>       | Standard Push button                                  |
| <b>Cartridge</b>           | Collect-type                                          |
| <b>Construction</b>        | Stainless steel exterior                              |
| <b>Drive</b>               | Air driven                                            |
| <b>Gear Ratio</b>          | 1:1                                                   |
| <b>LED Source Life</b>     | 1000 hours                                            |
| <b>Noise</b>               | 68 dB                                                 |
| <b>Rotation Speed</b>      | Up to 3,60,000 rpm                                    |
| <b>Torque Speed</b>        | 3,00,000 rpm                                          |
| <b>Tubing Connection</b>   | 6 holes                                               |
| <b>Type of Bearing</b>     | Ceramic                                               |
| <b>Weight</b>              | Approx. 180 g                                         |



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)