



High-Speed Dental Hand piece KDH-A201

www.nanomed.cc info@nanomed.cc

High-Speed Dental Hand piece KDH-A201

High-Speed Dental Hand piece KDH-A201 is a contra-angled air-driven dental drill instrument. High precision gears housed in miniature head design enables broader workspace. LED light illumination and single-point water spray with anti-retraction units improve scope of visibility.

Features:

Double optical LED
Anti-retraction valve mechanism
Enhanced grip control

Applications:

Used in dental clinics for preparation of tooth for various procedures primarily suitable for cavity preparation, removal of tooth structure, margin creation and finishing restorations.

Specifications :

Air Pressure	0.25 MPa ~ 0.27 MPa (4 hole) 0.2 MPa ~ 0.22 MPa (2 hole)
Applicable Bur Size	Diameter : 1.59 mm ~ 1.6 mm Length : 21 mm ~ 23 mm
Bur Insertion	Standard Push button
Cartridge	Collet-type
Construction	Stainless steel exterior
Drive	Air driven
Gear Ratio	1:1
Head Size	8.7 x 9.5 mm
LED Source Life	1000 hours
Noise	58 dB
Rotation Speed	3,80,000 rpm to 4,00,000 rpm
Spray	Single-point water spray
Tubing Connection	2 or 4 holes
Type of Bearing	Ceramic
Weight	Approx. 180 g



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc