



Diathermy Plate KDTP-A102

www.nanomed.cc info@nanomed.cc

Diathermy Plate KDTP-A102

KDTP-A102 is a bipolar neutral electrode plate with 130 cm² conductive area. Features secure connection of connector cable and electrode plate ensures user safety.

Features:

Neutral electrode plate
Type - Bipolar type (precise treatment)
Self -adhesivity ensures adherence
Lightweight and single use
Adult use

Applications:

Used in physical therapy, occupational therapy, rehabilitation department, surgical procedures, etc.

Specifications :

Conductive area	130 cm ²
Dimension	210 mm x 105 mm
Type	Bipolar
Use	Adult



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou,
Guangdong, China
Email: info@nanomed.cc | Website: www.nanomed.cc