



Infrared Ear Thermometer KETH-C100

www.nanomed.cc info@nanomed.cc

Infrared Ear Thermometer KETH-C100

KETH-C100 Infrared Ear Thermometer has a large screen display and ergonomic design that makes it simpler to use. Device has storage capacity up to 10 temperature reading.

Features:

Unlimited usage
Long life endurance
Fast response

Applications:

Used in clinics to measure the infrared energy radiated from ear canal.

Specifications :

Battery	2, AAA Battery
Dimension	134 x 70 x 34 mm
Measurement Time	1 s
Temperature Measuring accuracy	± 0.2 ?C
Temperature Measuring range	34 ?C ~ 44 ?C
Weight	75 gm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou,
Guangdong, China
Email: info@nanomed.cc | Website: www.nanomed.cc