



Dental X-ray KDX-B100

www.nanomed.cc info@nanomed.cc

Dental X-ray KDX-B100

Dental OPG X-ray KDX-B100 is a panoramic extra-oral imaging modality unit designed on the principle of curvature tomography. System utilizes oil immersed and self-cooling technology for longer sustenance. Micro-focus technology ensures clear image and accurate diagnosis.

Features:

Full 360° rotation
Micro focus technology
Automated Elevation motor
Maxillary sinus exam

Applications:

Used in dental clinics for dental X-ray imaging to reveal impacted wisdom teeth, maxillary arch, mandible, lower jaw or jawbone, maxillary sinuses, nasal sinuses and the whole 32 teeth area

Specifications :

Charger Input	AC : 220 V $\pm$ 10 %, 50 Hz ,1.8 kVA
Exposure Time	16 s
Film Size	150 $\times$ 300 mm
SSD (Source Skin Distance)	1120 mm
Tube Current	10 mA
Tube Focus	0.8 mm
Tube Voltage	90 kVp
Weight	228 kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc