



Otoscope KOO-A100

www.nanomed.cc info@nanomed.cc

Otoscope KOO-A100

Otoscope KOO-A100 is an optical handheld monocular illumination device for visual examination of ear-canal. Scratch-resistant high quality optical glass lens ensures clarity for detailed examination of anatomical structures. High intensity cool LED illumination bulb provides higher luminance with aperture selection promoting precise optical guidance.

Features:

- Focussed and homogeneous luminosity
- Precise magnification scope
- Scratch resistant glass lenses
- Washable & Re-usable Ear Tips
- Durable and comfortable handle

Applications:

Used in ENT, OPDs, examination clinics, audiology during rounds as well as for home-care use in pneumatic examination and magnification for detailed inspection of the in-ear anatomical structures and easy recognition during routine check-ups.

Specifications :

Batteries	Replaceable cells
Ear Specula sizes	2.4 mm 3 mm 4 mm 5 mm
Lens	Glass lens with 3 X magnification
Light source	LED bulb
Luminance	About 18.5 lumen
Patient groups	Infant to Adults



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc