



Fetal Doppler KFD-A101

www.nanomed.cc info@nanomed.cc

Fetal Doppler KFD-A101

Fetal Doppler KFD-A101 is a hand held ultrasound doppler with heart-rate display and audio output. Doppler probe ensures increased sensitivity with high fidelity. LCD display indicates fetal heart rate with graphs, volume and battery status. Built-in loudspeaker enables clear audibility of fetal heart sounds.

Features:

- High sensitivity doppler probe
- Low power consumption
- Compact battery powered operation
- Provision for earphones
- Backlit colour LCD display

Applications:

Used in gynecology clinics and obstetric department in hospitals to detect and monitor fetal heart rate for routine checkup.

Specifications :

Battery	Option A	Option B
	1.5 V AA battery	Ni-MH rechargeable battery
Dimension	135 × 95 × 35 mm	
Display	FHR and curve	
FHR accuracy	± 1 bpm	
FHR Measuring range	50 ~ 240 bpm	
FHR resolution	1 bpm	
Power consumption	1 W	
Screen	Color LCD	
Ultrasound frequency	2 MHz	
Ultrasound Intensity	10mW/cm ²	
Weight	500 g	



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc