



Knee Walker Scooter KNWS-A100

www.nanomed.cc info@nanomed.cc

Knee Walker Scooter KNWS-A100

KNWS-A100 is four wheeled scooter available in various color options built with adjustable height feature for handle and knee cushion. Steady and firm whole integral steel tube welded frame.

Features:

Steady and firm steel tube welded frame
Adjustable handle bar and knee cushion height
Foldable handle
Anti-slip back wheel band type brake
PU solid and wearable tire

Applications:

Used in hospitals and at homes. Suitable for patient with injury at lower limb.

Specifications :

Adjustable Handle Height range	31 to 41 inch
Dimension	780 × 420 × 775 to 1015 mm
Knee Pad Height range	17.5 to 22 inch
Material	Steel tubes coat with electroplates
Max Load	150 Kg
Weight	10.5 Kgs
Wheel Size	8 inch



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc