



Folding Power Wheelchair KFW-A100

www.nanomed.cc info@nanomed.cc

Folding Power Wheelchair KFW-A100

KFW-A100 is foldable and portable electrically-powered wheelchair easy to carry along and store for travel. Features detachable and washable seat cushion and backrest and flipping foot pedal. Wheelchair incorporates powerful motors responsible to easily operate wheelchair controller for user comfort.

Features:

Portable and foldable electric wheel chair

Dual motorized operation

Aluminum alloy frame

Detachable and washable seat cushion and backrest

Flipping foot pedal

Stable and durable black metal hub

Safe electromagnetic brake

Lithium battery and battery indicator

Applications:

Used in Hospitals and at Home for geriatric patient. Easily portable because of its folding nature.

Specifications :

Armrest height	220 mm
Back seat height	470 mm
Backrest height	360 mm
Battery	24 V / 10 Ah
Charging time	6 – 8 hour
Controller	24 V / 2 A
Dimension	870 × 610 × 920 mm
Driving distance	20 Km
Fold size	610 × 260 × 770 mm
Foot rest height	120 mm
Grade ability	8 – 12 °
Max load capacity	150 Kgs
Max speed	6 Km/hr
Motor	180 W brushless
Net weight	23 Kgs (with 1 battery)
Seat depth	400 mm
Seat height	520 mm
Seat width	400 mm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc