



Folding Power Wheelchair KFW-A101

www.nanomed.cc info@nanomed.cc

Folding Power Wheelchair KFW-A101

KFW-A101 is foldable and portable electrically-powered wheelchair easy to carry along and store for travel. Features detachable and washable seat cushion and backrest and flipping foot pedal. Wheelchair incorporates powerful motors responsible to easily operate wheelchair controller for user comfort.

Features:

- Portable and foldable electric wheel chair
- Dual motorized operation
- Aluminum alloy frame
- Detachable and washable seat cushion and backrest
- Flipping foot pedal
- Stable and durable black metal hub
- Safe electromagnetic brake
- Lithium battery and battery indicator

Applications:

Used in Hospitals and at Home for geriatric patient. Easily portable because of its folding nature.

Specifications :

Armrest height	250 mm
Back seat height	480 mm
Backrest height	400 mm
Battery	24 V / 6 Ah
Charging time	6 – 8 hour
Controller	24 V / 2 A
Dimension	960 × 580 × 930 mm
Driving distance	25 Km
Fold size	580 × 325 × 780 mm
Foot rest height	120 mm
Grade ability	8 – 12 °
Max load capacity	180 Kgs
Max speed	6 Km/hr
Motor	250 W brushless
Net weight	28 Kgs (with 2 battery)
Seat depth	430 mm
Seat height	460 mm
Seat width	480 mm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc