



Dental Chair Unit KDU-A201

www.nanomed.cc info@nanomed.cc

Dental Chair Unit KDU-A201

Dental Chair Unit KDU-A201 is a low-mounted smart dental unit equipped with advanced motors and fold-type seat for exact patient positioning and dentist stool. Inbuilt electromagnetic valve for smooth shift facility, wide-seat and soft leather upholstery ensures superior patient comfort. Anti-glare, medical grade LED light ensures precise illumination. Integrated cuspidor, 9 memory control system and instrument tray with high resolution LCD screen to set parameters for detailed set-up of chair and instruments.

Features:

- 9 Memory control program
- Folding type dental unit
- Wide seating with extra leg space
- Selectable suction apparatus
- Polyurethane upholstery (optional)
- Top-mounted unit (optional)

Applications:

Used in dental clinics for dental examination to reveal impacted teeth, perform common dental procedures like dental restorations, dental brace, extractions, periodontal surgery, filling and dental implantations.

Specifications :

Configuration	Low-mounted Folding
Dimension	1600 x 1060 x 1250 mm
Memory Programs	9 Control Programs
Power Voltage	220V-50Hz / 110V-60Hz
Upholstery	Leather
Weight	340 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou,
Guangdong, China
Email: info@nanomed.cc | Website: www.nanomed.cc