



## Operating Table KOT-A100

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

## **Operating Table KOT-A100**

Operating Table KOT-A100 features high density sponge in one piece with no seam, mattress is wrapped with artificial leather. It is incorporated with high quality motor without noise formation and foot treadle brake device for easy movement of operating table. Separate section for legs and it uses spring gas control for counterbalance (support)

### **Features:**

- Electrically controlled operation
- Adjustable height
- Trendelenburg and lateral tilt
- Base seat and upright column covers made up of stainless steel
- Panel control and remote control
- Head and leg section folding
- Back section folding
- Long life-span motor

### **Applications:**

Used in hospitals, operating room to carry out surgeries on patients.

## Specifications :

<b>Adjustable height (Min to Max)</b>	( 750 to 1000 ) $\pm$ 50 mm
<b>Back section fold downward</b>	? 20 °
<b>Back section fold upward</b>	? 45 °
<b>Head section folds downward</b>	? 90 °
<b>Head section folds upward</b>	? 40 °
<b>Input power</b>	300 VA
<b>Lateral tilt</b>	? 15 °
<b>Leg section folds downward</b>	? 90 °
<b>Leg section folds outward</b>	? 90 °
<b>Length</b>	2000 $\pm$ 50 mm
<b>Power supply</b>	220 V (AC) $\pm$ 10 % ; 50 Hz
<b>Reverse Trendelenburg</b>	? 35 °
<b>Trendelenburg</b>	? 25 °
<b>Waist lifting</b>	? 100 mm
<b>Width</b>	500 $\pm$ 20 mm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)