



Operating Table KOT-A101

www.nanomed.cc info@nanomed.cc

Operating Table KOT-A101

Operating Table KOT-A101 is a table top unit made up of double decker plastic plate used for surgical procedures. X-ray examination during operation is possible. Foot section can be positioned perpendicular to the leg section and is detachable. Raising, tilting and bending operation are controlled electrically.

Features:

- Adjustable height
- Freely raised, lowered and fixed within 250mm
- Trendelenburg and lateral tilt
- Head and leg section folding
- Back section folding

Applications:

Used in hospitals, operating room to carry out surgeries on patients.

Specifications :

Adjustable height (Min to Max)	750 to 1000 mm
Back section fold downward	? 15 °
Back section fold upward	? 45 °
Head section folds downward	? 90 °
Head section folds upward	? 45 °
Kidney bridge elevation	? 100 mm
Lateral tilt	? 15 °
Leg section opened	? 90 °
Length	2100 mm
Power supply	220 V (AC) ± 10 % ; 50 Hz
Reverse Trendelenburg	? 25 °
Trendelenburg	? 25 °
Width	480 mm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc