



Hydraulic Operating table KOT-A200

www.nanomed.cc info@nanomed.cc

Hydraulic Operating table KOT-A200

Operating Table KOT-A200 features high density sponge in one piece with no seam, mattress is wrapped with artificial leather. Controlled by hydraulic pneumatic systems to provide accurate positioning during major surgeries. Base of table is fixed and table can be rotated 360 degree. Carbon plastic plate on table surface allow easy X-ray examination.

Features:

Operating mode : Hydraulic system
Adjustable height, Trendelenburg and lateral tilt
High quality stainless steel rails on bed sides
Alloy aluminum fabricated table support
Base seat and upright column covering of stainless steel

Applications:

Used in hospitals, operating room to carry out surgeries on patients.

Specifications :

Adjustable height (Min to Max)	(740 to 1000) ± 50 mm
Head section folds downward	90 °
Head section folds upward	30 °
Horizontal rotary	360 °
Lateral tilt	25 °
Leg section folds downward	90 °
Leg section folds outward	180 °
Length	2000 ± 50 mm
Lumbar lower corner	120 °
Lumbar upper corner	200 °
Reverse Trendelenburg	20 °
Trendelenburg	20 °
Width	500 ± 20 mm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc