



## Multipurpose Operating Table KOT-B100

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

## **Multipurpose Operating Table KOT-B100**

KOT-B100 is multipurpose operating table designed such that all movements are manipulated from head. It provide comfort positions so that doctors can perform various surgeries and incorporates detachable foot section. Trendelenburg positions and left or right tilt are regulated by one single hand wheel. Oil pump controls the ascending and descending movement of the table.

### **Features:**

- Head control operating table
- Adjustable height
- Trendelenburg and lateral tilt
- Free height adjustment within 240 mm
- Allows 90° fold of foot section from leg section
- Hand wheel regulate back section and leg section movement
- Crank handle controls kidney bridge
- Leg section rotated by crank handle

### **Applications:**

Used in hospitals, operating room to carry out surgeries on patients.

## Specifications :

<b>Adjustable height (Min to Max)</b>	800 to 1040 mm
<b>Back section downward</b>	? 8 °
<b>Back section upward</b>	? 80 °
<b>Head section folds downward</b>	? 45 °
<b>Head section folds upward</b>	? 45 °
<b>Kidney section elevation</b>	? 120 mm
<b>Lateral tilt</b>	? 20 °
<b>Leg section folds downward</b>	? 90 °
<b>Leg section folds outward</b>	? 90 °
<b>Length</b>	2040 mm
<b>Reverse Trendelenburg</b>	? 50 °
<b>Trendelenburg</b>	? 30 °
<b>Width</b>	480 mm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)