



Hydraulic Operating Table KOT-C101

www.nanomed.cc info@nanomed.cc

Hydraulic Operating Table KOT-C101

Hydraulic Operating Table KOT-C101 is a multifunctional parturition table with hydraulic lift -n-lower mechanism of platform for safer operation. Modular table top sections allow quick set up and alignment. Streamlined cushioning with ergonomic design increases comfort.

Features:

- Anti-corrosive 304 stainless steel base
- Durable Y type sealing ring
- High density sponge mattress
- Rotary brake for ease of movement
- Extendable auxillary bed

Applications:

Used in hospitals, labour delivery room, operating room for performing gynecology operations, diagnosis, inspections and orthopaedic procedures on patients.

Specifications :

Adjustable Height (Min to Max)	(700-950) ± 50 mm
Auxillary Table	410 × 600 mm
Back Section	Upward: ? 50° Downward: ? 3°
Length	1900 ± 50 mm
Reverse Trendelenburg	20°
Seat Section	800 × 600 mm
Trendelenburg	10°
Width	600 ± 10 mm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc