



Halogen Surgical Light KOL-A200

www.nanomed.cc info@nanomed.cc

Halogen Surgical Light KOL-A200

KOL-A200 is an aperture operating light series featuring 12 bulbs with a power of 300 W which assists in daily surgical procedures. Wall mounted with minimum height of 2800 mm to 3200 mm for installation and uses halogen lamps as a light source for easy and clarified user observation.

Features:

- High illumination intensity
- Wall mounted lights
- Reduced shadow generation
- Eliminates overheating of surrounding

Applications:

Used in hospitals for surgeries.

Specifications :

Bulb power	300 W
Bulb voltage	24 V
Color temperature	4500 K ± 500 K
Illumination Intensity (at 1 m)	≥120,000 lux
Light source	Halogen bulb
Minimum height for installation	2800 to 3200 mm
Power input	350 VA
Power supply	AC 220 V±10 % 50 Hz



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China
Email: info@nanomed.cc | Website: www.nanomed.cc