



LED Surgical Light KOL-A400

www.nanomed.cc info@nanomed.cc

LED Surgical Light KOL-A400

KOL-A400 is shadow less operating light featuring LED cold light with ultrathin optical lenses suitable for surgical procedures. Digital step less regulated LED brightness and handle is adjustable to alter the light field size. Wall mounted with minimum height for installation of 2800 mm to 3200 mm.

Features:

- Light source : LED cold light
- Removable handle sheath
- Long life LED lamp
- Adapts switch power supply control voltage
- Spring arm for lamp head adjustment

Applications:

Used in hospitals as medical exam light, operating room lights.

Specifications :

Average bulb life	50000 hr.
Combined light depth of illumination	130 mm
Heat to light ratio	6
Illumination capacity	40000 lux ? 160000 lux / 40000 lux ? 160000 lux
Light field diameter	100 mm to 280 mm
Light head	550 mm
Operating distance	700 ? 1400 mm
Optional Color temperature	(3500 / 4000 / 4500) $\pm$ 500K
Power input	180 VA
Power supply	100 / 240 V 50 / 60 Hz



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc