



Connecting Stretcher KCS-A101

www.nanomed.cc info@nanomed.cc

Connecting Stretcher KCS-A101

KCS-A101 is emergency connecting trolley for operating room. Stretchers can be moved for easy transport of patient from outside to inside the operating room when two trolleys are at the same level. Stretcher can be locked automatically with the trolley. Two trolleys with rise and fall guide wheel are joined together but up and down separately. Backrest function is controlled by spring.

Features:

- Aluminum alloy trolley structure
- Aluminum alloy railway to connects 2 stretchers
- Two steel flat side rails
- Steel bottom platform
- Central locking function
- IV pole, PU (polyurethane) cover water proof mattress

Applications:

Used in hospital operating rooms to transfer patient from table to anesthesia.

Specifications :

Backrest upward angle	75 °
Bottom platform	Steel material
Caster diameter	150 mm
Dimension	3650 × 640 × 640-970 mm
Mattress thickness	80 mm
Side rails	Steel material
Variable height	640 to 970mm
Weight capacity	150 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc