



## Electric Stretcher KES-A100

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

## Electric Stretcher KES-A100

KES-A100 is a functional stretcher with a provision of five movements namely backrest and footrest movement, height adjustment, trendelenburg and anti-trendelenburg position control. Equipped with electric motors for effortless lifting operation and electric control panel for operational monitoring and controlling.

### Features:

- Electric motors
- Spring controlled back section, manually controlled foot rest
- Aluminum alloy handrails provide safety protection
- CPR handle at head side for quickly lowering the head section
- Steel bottom platform
- Entire bed platform with x-ray cassette for imaging
- Electric control pedal
- Rise and fall fifth guide wheels
- PU cover water proof mattress prevent dust mites and allergies
- IV pole, IV pole holders
- Central locking function

### Applications:

Used in hospitals, nursing homes for long stay.

## Specifications :

|                              |                           |
|------------------------------|---------------------------|
| <b>Anti-trendelenburg</b>    | 15 °                      |
| <b>Backrest upward angle</b> | 75 °                      |
| <b>Dimension</b>             | 2100 × 850 × 720 – 950 mm |
| <b>Footrest upward angle</b> | 45 °                      |
| <b>Mattress thickness</b>    | 80 mm                     |
| <b>Trendelenburg</b>         | 15 °                      |
| <b>Variable height</b>       | 720 – 950 mm              |
| <b>Weight capacity</b>       | 250 Kgs                   |



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)