



Electric Stretcher KES-A102

www.nanomed.cc info@nanomed.cc

Electric Stretcher KES-A102

KES-A102 is a functional stretcher with a provision of two movements namely backrest and height adjustment. Equipped with electric motors for effortless lifting operation and electric control panel for operational monitoring and controlling.

Features:

- Electric motors
- Spring controlled back section
- ABS handrails provide safety protection
- CPR handle at head side for quickly lowering the head section
- PP (polypropylene) bottom platform
- Steel frame
- Entire bed platform with x-ray cassette for imaging
- Electric control pedal
- Rise and fall fifth guide wheels
- IV pole
- Central locking function

Applications:

Used in hospitals, nursing homes for long stay.

Specifications :

Backrest upward angle	75 °
Dimension	1930 × 640 × 540 – 840 mm
Mattress thickness	60 mm
Variable height	540 – 840 mm
Weight capacity	150 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou,
Guangdong, China
Email: info@nanomed.cc | Website: www.nanomed.cc