



Hydraulic Stretcher KHS-A102

www.nanomed.cc info@nanomed.cc

Hydraulic Stretcher KHS-A102

KHS-A102 is a functional stretcher with a provision of five movements namely backrest and footrest movement, height adjustment, trendelenburg and anti-trendelenburg position control. Equipped with 2 separate hydraulic pumps for lifting operation and Hydraulic control panel for monitoring and controlling operation. Entire bed board features x-ray cassette for imaging with a monitor tray at foot end.

Features:

- Two separate sets of hydraulic pumps
- Spring controlled back section, manually controlled foot rest
- Aluminum alloy handrails provide safety protection
- CPR handle at head side for quickly lowering the head section
- Entire bed board features x-ray cassette for imaging
- Cold-rolled steel tube frames
- Steel, electrostatic powder coated mattress surface
- Rise and fall fifth guide wheels ensure shock proof
- PU cover water proof mattress with straps prevent dust mites and allergies
- IV pole, IV pole holders
- ABS base cover with oxygen cylinder holder
- Central locking function

Applications:

Used in hospitals and clinics to move the patient's and to lift the patient with spinal cord injury.

Specifications :

Anti-trendelenburg	12 °
Backrest upward angle	75 °
Castor diameter	200 mm
Dimensions	2100 mm × 850 mm × 630 - 950 mm
Footrest upward angle	45 °
Mattress thickness	60 mm
Trendelenburg	12 °
Variable height	630 to 950 mm
Weight capacity	250 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc