



Manual Transfer Stretcher KMS-A100

www.nanomed.cc info@nanomed.cc

Manual Transfer Stretcher KMS-A100

KMS-A100 functional stretcher with a provision of two movements namely backrest and height adjustment. It allows manual transfer of patient onto another location. Back section movement is controlled by spring. Fifth rise and fall guide wheel features shock proof system.

Features:

- Cold-rolled steel tube frame
- PP (Polypropylene) bottom platform
- CPR handle at head side for quickly lowering the head section
- Rise and fall guide wheels
- IV pole, IV pole holders
- ABS bed board with oxygen cylinder holder
- Water proof mattress with straps

Applications:

Used in hospital, nursing home and ambulances.

Specifications :

Back rest upward angle	75 °
Castor diameter	150 mm
Dimension	1930 × 640 × 500 - 900 mm
Mattress	60 mm thick
Transfer operation	Manual
Variable height	500 to 900 mm
Weight capacity	250 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc