



Manual Transfer Stretcher KMS-A101

www.nanomed.cc info@nanomed.cc

Manual Transfer Stretcher KMS-A101

KMS-A101 functional stretcher with a provision of height adjustment. It allow manual transfer of patient onto another location. Features four wheels and artificial leather mattress.

Features:

Stainless steel frame
IV pole
Four wheels
Artificial leather mattress

Applications:

Used in hospital, nursing home and ambulances.

Specifications :

Dimension	1900 × 630 × 600 - 800 mm
Mattress	Artificial leather
Transfer operation	Manual
Variable height	600 to 800 mm



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou,
Guangdong, China
Email: info@nanomed.cc | Website: www.nanomed.cc