



## **CPAP Machine KCM-A100**

[www.nanomed.cc](http://www.nanomed.cc) [info@nanomed.cc](mailto:info@nanomed.cc)

## **CPAP Machine KCM-A100**

CPAP Machine KCM-A100 is a noninvasive ventilation system equipped with LCD screen for display of patient trends, CPAP pressure. Uniform humidification ensures optimal respiratory comfort. Patient respiration automatically triggers work mode and expiratory pressure relief (EPR) function traces every breath for safe and reliable performance. Internal memory records case details for data transmission via USB port. Integrated automatic leakage and altitude compensation for improved accuracy.

### **Features:**

- 3.5 inch TFT LCD screen
- 128x64 LCD resolution
- Accurate and reliable performance
- Single touch operation
- Quiet operation
- Modes: CPAP and CPAP+EPR
- Patient data recording
- Timely indicative visual alarm

### **Applications:**

Used in clinical or primary healthcare facilities to treat or alleviate respiratory ailments such as obstructive sleep apnea syndrome, chronic obstructive pulmonary disease, acute lung injury, pneumonia, ventilator weaning, etc.

## Specifications :

<b>Adjustable Uniform Humidification System</b>	Yes
<b>Automatic Altitude Adjustable</b>	Yes
<b>Automatic Leakage Compensation</b>	Yes
<b>Automatic ON/OFF</b>	Yes
<b>Dimensions</b>	255 × 170 × 112 mm
<b>Display</b>	3.5'' TFT LCD
<b>EPR</b>	2-4 cm H2O
<b>Input</b>	110 V - 240 V AC, 50 Hz/60 Hz
<b>Mode</b>	CPAP, CPAP + EPR
<b>Noise Level</b>	30 dBA
<b>Pressure</b>	4-20 cm H2O
<b>Ramp Duration</b>	0-60 min
<b>Replaceable Water Chamber</b>	Yes
<b>Resolution</b>	128 X 64 LCD
<b>USB Connectivity</b>	Yes
<b>Weight</b>	1.8 Kg



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: [info@nanomed.cc](mailto:info@nanomed.cc) | Website: [www.nanomed.cc](http://www.nanomed.cc)