



Oxygen Concentrator KOC-A101

www.nanomed.cc info@nanomed.cc

Oxygen Concentrator KOC-A101

Oxygen Concentrator KOC-A101 is a 5L continuous flow portable medical oxygen concentration unit. Equipped with intelligent oxygen delivery system and clinically advanced medical grade molecular sieve with micro-particle entrapment filter ensures delivery of pure oxygen for safe inhalation as per oxygen flowrate control. Operates on battery as well as power supply to provide hours of consistent oxygen supply with minimal noise at low maintenance requirements.

Features:

- Capacity: 5L
- Touch screen control
- Silent operation
- High concentration of oxygen
- Variable flow-rate settings

Applications:

Used in primary healthcare facilities, ambulatory oxygen therapy as well as for point of care treatments at prescription for administrating point of care delivery of supplemental oxygen supply to maintain targeted SpO2 levels consistently all day to patients with respiratory disorders.

Specifications :

Capacity	5 L
Dimensions	335 x 215 x 375 mm
Flow Range	1 L/min – 5 L/min
Noise	? 60 dB
Operation	Continuous Flow
Operational temperature	5 °C to 40 °C
Oxygen Concentration	30 % – 90 %
Power Output	? 200 W
Power supply	220 V/110 V, 50 Hz/ 60 Hz
Weight	6.5 Kgs



Guangzhou Nano Medical Equipment Co., Ltd. 407, No. 1 Qingbu Street, Xinya Street, Huadu District, Guangzhou, Guangdong, China

Email: info@nanomed.cc | Website: www.nanomed.cc